EN	Instructions for Use
DE	Gebrauchsanweisung
FR	Mode d'emploi
ES	Instrucciones de uso
IT	Istruzioni per l'uso
SV	Bruksanvisning
NO	Bruksanvisning
DA	Brugsanvisning
FI	Käyttöohjeet
EL	Οδηγίες χρήσης
NL	Gebruiksaanwijzing
PL	Instrukcja użytkowania
PT	Instruções de Utilização
RU	Руководство по эксплуатации
EE	Kasutusjuhend
SL	Navodila za uporabo
CZ	Návod k použití
TR	Kullanma talimatları
HR	Upute za uporabu
HU	Használati útmutató
SR	Uputstvo za upotrebu
RO	Instrucțiuni de utilizare
AR	تعليمات الاستخدام
KO	사용 설명서
JP	取扱説明書
ZH	使用说明
BG	Инструкции за употреба

sedaconda®ACD

Anaesthetic Conserving Device

Inspiratory limb usage

INDEX //

EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Mode d'emploi	8
ES	Instrucciones de uso	10
IT	Istruzioni per l'uso	12
SV	Bruksanvisning	14
NO	Bruksanvisning	16
DA	Brugsanvisning	18
FI	Käyttöohjeet	20
EL	Οδηγίες χρήσης	22
NL	Gebruiksaanwijzing	24
PL	Instrukcja użytkowania	26
PT	Instruções de Utilização	28
RU	Руководство по эксплуатации	30
ET	Kasutusjuhend	32
SL	Navodila za uporabo	34
CS	Návod k použití	36
TR	Kullanma talimatları	38
HR	Upute za uporabu	40
HU	Használati útmutató	42
SR	Uputstvo za upotrebu	44
RO	Instrucțiuni de utilizare	46
AR	تعليمات الاستخدام	48
KO	사용 설명서	50
JP	取扱説明書	52
ZH	使用说明	54
BG	Инструкции за употреба	56

1. INTENDED USE

This Instructions For Use (IFU) describes the use of Sedaconda® ACD for Small Tidal Volumes and the connection of the Sedaconda® ACD to the inspiratory port of the ventilator. This placement is intended for and enables the delivery of inhaled anaesthetics to patients with tidal volumes of 30-200 ml. The advantage of the inspiratory side placement is that no dead space is added to the breathing circuit.

In contrast to the standard placement of the Sedaconda® ACD (between the ventilator/Y-piece and the patient) the inspiratory side placement only utilises the evaporator function of the Sedaconda® ACD and there is no reflection of inhaled anaesthetic. Therefore, higher inhaled anaesthetic pump rates can be expected than with standard placement, despite lower tidal volumes.

Administration of isoflurane and sevoflurane using Sedaconda® ACD should only be done in a setting fully equipped for the monitoring and support of respiratory and cardiovascular function and by persons specifically trained in the use of inhalational anaesthetic drugs and the recognition and management of the expected adverse effects of such drugs, including respiratory and cardiac resuscitation. Such training must include the establishment and maintenance of a patient airway and assisted ventilation.

For more information regarding standard placement, see Sedaconda® ACD IFU.

PRINCIPLES OF OPERATION

The Sedaconda® ACD consists of a plastic housing with an agent line for the continuous delivery of isoflurane or sevoflurane from a syringe pump to the miniature vaporizer where any clinical dosage is immediately vaporized.

2. IMPORTANT USER INFORMATION

2.1 Carefully read these instructions before using Sedaconda® ACD and note the following

GENERAL WARNINGS

- Do Not use an Sedaconda® ACD if the integrity of the package is breached or if packaging is visibly damaged.
- Do Not Use an additional standard Heat and Moisture Exchanger (HME) filter when the Sedaconda® ACD is used in the inspiratory side placement. Use of an HME filter will increase dead space in the circuit and add resistance due to accumulation of water. Instead active humidification of the respiratory gas is mandatory. Without active humidification the patient will be subject to dry medical gas without humidification.
- Always place the active humidification device below the Sedaconda® ACD, to avoid accumulation of condensate, with the black face uppermost.
- Do Not re-connect a used Sedaconda® ACD that has been disconnected and unattended for any reason for any length of time. Always use a new one. There is a risk of losing control over the concentration of volatile agent in the Sedaconda® ACD, specifically a risk of over dosing by unintentional filling from the syringe.
- Do Not use the Sedaconda® ACD gas sampling port, since it is not giving correct gas-measurements in the inspiratory side placement.
- Always stop the syringe pump if disconnecting the Sedaconda® ACD.
- Do Not use the bolus or purge function on the syringe pump unless programmed according to hospital protocol.
- Sedana recommends the use of pre-programmable function on the syringe pump when delivering a bolus in order to minimise the risk of overdosing.
- Do Not fold or clamp the agent line.
- Do Not seal the connector on the ventilator side except at disposal of Sedaconda® ACD
- Do Not use Sedaconda® ACD with jet or oscillation ventilation.
- Re-processing of medical devices intended for single use only may result in degraded performance or a loss of functionality e.g. resistance to breathing might increase. This product is not designed to be cleaned, disinfected, or sterilized.
- Only use CE-approved ventilators complying with applicable requirements, including standard ISO 80601-2-12. Sedaconda® ACD can be used with all conventional ventilator modes except oscillator mode for intubated patients. Use ventilator circuits compatible with anaesthetic agents. Use only with ventilators with an accessible exhaust.
- Only use CE-approved syringe pumps complying with applicable requirements, including standard ISO 60601-2-24. The pump must be programmable for Becton Dickinson Plastipak/Sherwood Monoject 50 or 60 ml syringes. The highest switch-off pressure should be applied, to eliminate pressure alarms due to the narrow lumen of the anaesthetic agent line.
- Anaesthetic gases should be monitored with a CE-approved gas analyser, which complies with its applicable requirements and with the specifications of standard ISO 80601-2-55.
- Sedaconda® ACD is approved for the delivery of volatile agents (VA). Sedana specifically recommends the use of isoflurane or sevoflurane, not desflurane or other volatile agents.
- Use polycarbonate-based components with precaution. Components may become degraded or undergo stress cracking, if used in the patient breathing circuit in the presence of the anaesthetic gases isoflurane or sevoflurane.


Symbol	Description
WARNING!	Indicates a condition which if not followed exactly may cause harm to patient or a user. Do not proceed until the instructions are clearly understood and all stated conditions are met.
	Indicates a condition, which if not followed exactly may cause harm to the product or equipment. Do not proceed until the instructions are clearly understood and all stated conditions are met.
NOTE!	Indicates information important for optimal use of the product.
	Single Patient Use
	Read the Instruction for Use carefully before use
	Not for IV use
	Ventilator
V_T	Tidal Volume
	Change Every 24 Hours
	Patient / Lungs
	Gas Measuring Device
	Keep stored away from direct light sources

3. PARTS DIAGRAM

The materials needed for the assembly (fig 1)

- Syringe pump and Sedaconda® Syringe
- Anaesthetic gas analyser
- Gas sampling line/Nafion line
- Airway connector
- Ventilator
- Sedaconda® ACD-L or Sedaconda® ACD-S
- Gas scavenging system
- Active humidifier

Fig. 1

4. SETUP

4.1. Connecting scavenging

Residual anaesthetic agent scavenging is recommended when using the Sedaconda® ACD. Connect active or passive scavenging according to the manufacturer's instructions for use. During inspiratory side placement additional moisture accumulates in the filter, which causes it to fill quicker. The FlurAbsorb should therefore be changed after 5 syringes (à 50 ml), 24 hours or on demand when resistance increase. For more information see the IFU for FlurAbsorb.

4.2. Connecting the anaesthetic gas analyser

- If the gas analyser requires a water trap or a special connection/adaptor, use a new connector (check the replacement intervals in the IFU of the gas analyser).
- Connect the gas analyser exhaust to the scavenging system.
- Switch on the gas analyser.
- Connect the airway connector between the Y-piece and the endotracheal tube.
- Connect the Nafion line to the airway connector.
- Connect the gas sampling line between the Nafion line and the gas analyser.


4.3. Filling the syringe

- Always store isoflurane and sevoflurane at room temperature.
- Do Not use an Sedaconda® Syringe that has been pre-filled and stored for longer than 5 days. Sedaconda® Syringe is for single patient use.
- Do Not fill the Sedaconda® Syringe without the Filling Adaptor. Attempting to fill the syringe without the adaptor or with a different adaptor could lead to accidentally using the wrong medication with Sedaconda® ACD.
- Open the bottle and screw the Filling Adaptor on to the bottle.
- Unscrew the red cap from the red top of the Sedaconda® Syringe.
- Fill the Sedaconda® Syringe with 10 to 20 ml of air.
- Connect the Sedaconda® Syringe tightly to the Filling Adaptor.
- Turn the bottle upside down. Fill the Sedaconda® Syringe slowly to avoid excessive bubbles.
- Note: To avoid bubbles, try slowly moving the Syringe plunger back and forth.
- Unscrew the Syringe from the filling adaptor.
- Remove any air from the Sedaconda® Syringe and close with the red cap.
- Label the Sedaconda® Syringe with information regarding anaesthetic agent and date of filling or label according to Hospital protocol.

4.4. Connecting the Sedaconda® ACD

- Always use active humidification when the Sedaconda® ACD is placed on the inspiratory side.
- Always place the Sedaconda® ACD above the active humidifier device.
- Set up the active humidifier circuit according to the manufacturer's instructions.

- Remove the red cap from the Sedaconda® ACD and the purple label from the gas sampling port and close it.
- Connect the Sedaconda® ACD to the inspiratory port of the ventilator.
- Position the Sedaconda® ACD with the black side facing up and the patient side connection point sloped down towards the active humidifier, with an angle of about 45 degrees. If not possible use a short flex tube.
- Push the flexible extension tube on to the patient side connection point on the Sedaconda® ACD, and the other side of the tube to the humidifier and subsequently to the patient.

5. SETUP

5.1. Priming and start of therapy

- ▶ Check that all connections are secure before starting therapy
- ▶ Do Not kink or clamp the anaesthetic agent line. This might damage the line. The line contains a valve that prevents medication from escaping.
- ▶ Do Not manually prime or give manual bolus doses, instead always use the syringe pump.
Manual priming and blousing may result in an overdose.
- Place the filled Sedaconda® Syringe in the syringe pump.
- Open the red syringe cap and connect the Sedaconda® ACD anaesthetic agent line to the Sedaconda® syringe.
- Program the syringe pump to prime the anaesthetic agent line with a bolus of 1.2 ml. Never prime manually.
- Start the syringe pump at 2 ml/hour and titrate the clinical dosage according to patient requirements.

5.2. Dose titration

- Titrate the dose according to clinical evaluation. The use of and dose of concomitant medications should be reviewed when starting treatment with Sedaconda® ACD.
- Increase the syringe pump rate in steps of 0.5 to 1 ml/hour until the desired FET% has been reached.
- Monitor the patient clinical status and hemodynamics closely to determine whether bolusing or dose changes are needed.
- If necessary, program the syringe pump to give small bolus doses of 0.2 to 0.3 ml. Never use the purge function or manual bolus. Caution: the age and size as well as the condition of the patient need to be taken into account.

5.3. Ending therapy

There are two options to end the therapy.

For rapid wake-up, stop the syringe pump. The FET% directly starts to drop – the patient will normally wake up within minutes unless other drugs are given.

For slow wake-up, reduce the syringe pump rate stepwise according to clinical needs. The FET% will be reduced according to the reduced pump rate.

6. REPLACING AN EMPTY SEDACONDA® SYRINGE AND SEDACONDA® ACD DEVICE

In order to reduce the risk of inadvertent awakening, the Sedaconda® ACD and Sedaconda® Syringe should be changed without delay. Drug concentration will drop as soon as the syringe pump is stopped.

6.1. Changing the Sedaconda® Syringe

- ▶ Change the Sedaconda® Syringe after one use, as soon as it empties.
- Prepare a new Sedaconda® Syringe according to the steps above, part 4.3.
- Stop the syringe pump.
- Disconnect the anaesthetic agent line from the used Sedaconda® Syringe and close the Sedaconda® Syringe with the red cap.
- Remove the used Sedaconda® Syringe from the syringe pump.
- Place the filled Sedaconda® Syringe in the syringe pump and unscrew the red cap.
- Connect the Sedaconda® ACD anaesthetic agent line to the new Sedaconda® Syringe.
- Start the syringe pump with the same rate as before.
- Dispose of the used Sedaconda® Syringe according to hospital protocols.

6.2. Changing the Sedaconda® ACD

- ▶ Change the Sedaconda® ACD after 24 hours of use, or earlier in the event of abnormal secretions or blockages. The Sedaconda® ACD is for single patient use only.
- Prepare a new Sedaconda® ACD.
- Stop the syringe pump.
- Disconnect the anaesthetic agent line from the Sedaconda® Syringe and close the Sedaconda® Syringe with the red cap.
- Disconnect the used Sedaconda® ACD from the ventilator inspiratory port side first.
- Disconnect the used Sedaconda® ACD from the flexible extension tube side next.
- Insert the new Sedaconda® ACD by connecting at the inspiratory port first and then to the flexible extension tube.
- Connect the anaesthetic agent line to the Sedaconda® Syringe in the syringe pump.
- Program the syringe pump to prime the anaesthetic agent line with a bolus of 1.2 ml. Never prime manually.
- Start the syringe pump with the same rate as before.

7. DISPOSAL OF PRODUCTS

Syringes with remaining liquid anaesthetic must be disposed of according to local routines for special waste. All other products can be disposed of in the standard hospital waste.

Recommended replacement intervals:

Item	Changing interval
Sedaconda® ACD-L Sedaconda® ACD-S	Single patient use Changed every 24 hour or earlier if needed
Sedaconda® syringe	Single use
Nafion line Gas sampling line	Single patient use, with a maximum length use of 7 days Replaced as needed in accordance to hospital hygiene regulations
FlurAbsorb	Multi patients use Capacity of up to 5 syringes (à 50ml), 24 hours use or when signs of increased resistance are detected
FlurAbsorb Accessory Kit	Single patient use
Single Use Filling Adapter Isoflurane/Sevoflurane Standard screw top bottles	Single use
Water trap	Single patient use Max days of use according to the manufacturer's instructions

8. PROCEDURES AFFECTING DRUG DELIVERY

Procedures that may affect drug delivery include:

- Endotracheal suctioning
- Bronchoscopy
- Nebulisation

For endotracheal suction, there are two possible ways, either with closed suction system or adapter with bronchoscopy cap. If clinically needed a bolus of inhaled anaesthetic can be given shortly before suctioning.

During bronchoscopy minor leakage may occur. Open only the smallest opening of the bronchoscopy adapter in order to minimise leakage. Inhaled anaesthetic delivery is reduced by the bronchoscope. Additional medications should be considered to facilitate the procedure.

Ambient concentration measurements during these procedures have shown transient elevations but these elevations are short-lived and below recommended exposure limits, indicating minimal risks for staff.

During nebulisation, the nebuliser should be placed in the breathing circuit according to the manufacturer's instructions for use. Use the device according to normal daily procedures.

WARNING! Repeated nebulisations may increase the flow resistance of the Sedaconda® ACD. Pay attention to signs of occlusions Sticky nebulised drugs (e.g. acetylcysteine, colistin and amphotericin B) may increase the resistance of the filter and require early replacement of the Sedaconda® ACD.

TECHNICAL SPECIFICATION	100 ML	50 ML
Anaesthetic Agents	Only use room temperature sevoflurane (18°-25°C) and isoflurane (18°-30°C)	
Syringe	Only use the Sedaconda® syringe	
Stability of filled syringes	5 days	
Tidal volume working range	350-1200 ml at Y-piece & ET-tube	200-800 ml at Y-piece & ET-tube
	30-200 ml on inspiratory placement	30-200 ml on inspiratory placement
Sedaconda® ACD dead space	Approx. 100 ml	Approx. 50 ml
Resistance to gas flow at 60 l/min	2.5 cm H ₂ O (250 pa)	3.0 cm H ₂ O (300 pa)
Moisture loss	5 mg/l (@ 0.75L X 12 bpm)	5 mg/l (@ 0.5L X 15 bpm)
	7 mg/l (@ 1.0L x 10 bpm)	6 mg/l (@ 0.75L x15 bpm)
Filter capacity:	Bacterial filtration Viral Filtration	99,867 % 99,76 %
Weight	50 g	
Agent Line Length	2.2 m	
Connectors (According to ISO 5356)	15F/22M-15M	
Gas Sampling Port	Female Luer Lock	

In case of a serious incident please contact the competent authority and the legal manufacturer Sedana Medical Ltd.

 Sedana Medical Ltd., Unit 2A The Village Centre,
Two Mile House, Naas, Co. Kildare, W91 PWH5, Ireland

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