EN	Instructions for Use
DE	Gebrauchsanweisung
FR	Mode d'emploi
ES	Instrucciones de uso
IT	Istruzioni per l'uso
SV	Bruksanvisning
NO	Bruksanvisning
DA	Brugsanvisning
FI	Käyttöohjeet
EL	Οδηγίες χρήσης
NL	Gebruiksaanwijzing
PL	Instrukcja użytkowania
PT	Instruções de Uso
RU	Руководство по эксплуатации
ET	Kasutusjuhend
SK	Návod na použitie
SL	Navodila za uporabo
CS	Návod k použití
TR	Kullanma talimatları
HR	Upute za uporabu
HU	Használati útmutató
SR	Instrukcije za upotrebu
RO	Instrucțiuni de utilizare
AR	تعليمات الاستخدام
KO	사용 설명서
JP	取扱説明書
ZH	使用说明
BG	Инструкции за употреба

Sedaconda® ACD-S  
Sedaconda® ACD-L  
Sedaconda® Syringe

Administration system for  
anaesthetic agents

## INDEX //

EN	Instructions for Use	4	EN
DE	Gebrauchsanweisung	6	DE
FR	Mode d'emploi	8	FR
ES	Instrucciones de uso	10	ES
IT	Istruzioni per l'uso	12	IT
SV	Bruksanvisning	14	SV
NO	Bruksanvisning	16	NO
DA	Brugsanvisning	18	DA
FI	Käyttöohjeet	20	FI
EL	Οδηγίες χρήσης	22	EL
NL	Gebruiksaanwijzing	24	NL
PL	Instrukcja użytkowania	26	PL
PT	Instruções de Uso	28	PT
RU	Руководство по эксплуатации	30	RU
ET	Kasutusjuhend	32	ET
SK	Návod na použitie	34	SK
SL	Navodila za uporabo	36	SL
CS	Návod k použití	38	CS
TR	Kullanma talimatları	40	TR
HR	Upute za uporabu	42	HR
HU	Használati útmutató	44	HU
SR	Instrukcije za upotrebu	46	SR
RO	Instrucțiuni de utilizare	48	RO
AR	مادخت سالا تامي لعنت	50	AR
KO	사용 설명서	52	KO
JP	取扱説明書	54	JP
ZH	使用说明	56	ZH
BG	Инструкции за употреба	58	BG

1. INTENDED USE

Sedaconda ACD (Anaesthetic Conserving Device) is intended for administering isoflurane or sevoflurane to invasively ventilated patients. Administration of isoflurane and sevoflurane using Sedaconda ACD should only be done in a setting fully equipped for the monitoring and support of respiratory and cardiovascular function and by persons specifically trained in the use of inhalational anaesthetic drugs and the recognition and management of the expected adverse effects of such drugs, including respiratory and cardiac resuscitation. Such training must include the establishment and maintenance of a patient airway and assisted ventilation. The Sedaconda ACD is intended for single-use only and needs to be replaced every 24 hours or when needed e.g. at unexpected events such as sudden blockage of the airways because of secretion etc.


The instructions contained in this document apply to Sedaconda ACD-S, Sedaconda ACD-L and Sedaconda Syringe.

2. PRINCIPLES OF OPERATION

The Sedaconda ACD consists of a plastic housing with an agent line for the continuous delivery of isoflurane or sevoflurane from a syringe pump to the miniature vaporizer where any clinical dosage is immediately vaporized. During continued breathing the volatile anaesthetic agent is captured and reflected by the activated carbon filter.

The Sedaconda ACD is available in two sizes: Sedaconda ACD-L (100ml dead space) and Sedaconda ACD-S (50ml dead space). The dead space of 100 ml or 50 ml needs to be considered for all patients and CO<sub>2</sub> needs to be carefully monitored. Adjustments to CO<sub>2</sub> can be achieved by optimising the ventilator parameters. In addition, the Sedaconda ACD is an effective heat and moisture exchanger and it includes an efficient bacterial/viral filter.

Cross Section Drawing of the Sedaconda ACD

3. IMPORTANT USER INFORMATION

3.1 Carefully read these instructions before using Sedaconda ACD and note the following

GENERAL WARNINGS

- Do not use desflurane.
- Do not re-connect a used Sedaconda ACD that has been disconnected and unattended for any reason for any length of time. Always use a new one.
- Only use and store isoflurane at 18-30°C.
- Only use and store sevoflurane at 18-25°C.
- Only use anaesthetic agent with the same temperature as the treatment room temperature.
- Do not use a Sedaconda ACD if the integrity of the package is breached or if packaging is visibly damaged.
- Always stop the syringe pump if disconnecting the Sedaconda ACD.
- Do not prime the agent line manually. Always use the syringe pump.
- Position the patient side connector of the Sedaconda ACD lower than the ventilator side, to avoid accumulation of condensate, the arrows on the device indicate the correct orientation.
- Do not use the bolus or flush function on the syringe pump unless programmed according to hospital protocol.
- Do not fold or clamp the agent line.
- Do not use Sedaconda ACD with jet or oscillation ventilation.
- Do not use active humidification together with Sedaconda ACD in the standard setup/ placement. See section 11 for alternative placement setup.
- Do not use Sedaconda ACD on patients with copious secretions in the standard setup/ placement.
- Re-processing of medical devices intended for single-use only may result in degraded performance or a loss of functionality e.g. resistance to breathing might increase. This product is not designed to be cleaned, disinfected, or sterilized.
- Never seal the connector on the ventilator side except at disposal of Sedaconda ACD.
- Only use isoflurane and sevoflurane according to their respective SmPC documents.
- Always consider the dead space of the device vs tidal volume when ventilating the patient and carefully monitor CO<sub>2</sub> level.
- Use the gas monitor port according to section 4.3 "Patient monitoring". Connection to other devices may result in patient hazard.

Symbol	Description				
	Not for IV use			<b>MD</b>	Medical Device
	Ventilator	$V_T$	Tidal Volume		Change Every 24 Hours
	Patient / Lungs		Gas Measuring Device		Keep stored away from direct light sources
	PVC free		phthalates free		Natural rubber latex free

4. ADDITIONAL EQUIPMENT REQUIRED (FIG 1)

Only medical devices which bear the CE mark and which comply with its applicable international standards, may be used. The following equipment should be available when using the Sedaconda ACD and Sedaconda Syringe:

- Syringe pump with settings for BD Plastipak or Monoject Sherwood 50/60 ml syringes

- Anaesthetic gas monitor, which displays concentrations of CO<sub>2</sub> and anaesthetic gases
- Ventilator
- Applicable adapters for isoflurane and sevoflurane use
- Gas scavenging system

4.1 Sedaconda Syringe

The Sedaconda Syringe is validated under the BD Plastipak 50ml and Monoject Sherwood 50ml settings. The Sedaconda Syringe contains a non-standard, non-luer coupling designed to fit the Sedaconda ACD agent line connector. The syringe can be pre-filled and stored up to 5 days in a dark environment at room temperature. Make sure that the syringe is safely closed.

4.2 Syringe pumps

Use only CE-labeled syringe pumps, which comply with its applicable requirements, in particular with the specifications of standard EN 60601-2-24, and which are programmable pumps with settings for Becton Dickinson Plastipak or Sherwood Monoject 50/60 ml syringes. Ensure that the syringe pump is at or below the level of the Sedaconda ACD.

4.3 Patient monitoring

During treatment with Sedaconda ACD a gas analyser, which displays Fet (end-expiratory concentration) representing the alveolar concentration of volatile anaesthetics, should be available for measurements of gas concentration. Measurement of Fet is recommended upon starting therapy and informative when using neuromuscular blocking agents. The Fi concentration should not be used.

Both side stream and main-stream gas analysers can be used with Sedaconda ACD. The anaesthetic gas monitor should be CE-labeled complying with its applicable requirements, in particular with the specifications of standard EN ISO 80601-2-55.

4.4 Ventilator

Use only CE-labeled ventilators which comply with its applicable requirements, in particular with the specifications of standard ISO 80601-2-12. Sedaconda ACD can be used on all conventional modes but not on oscillator mode for intubated patients. Use ventilator circuits suitable for use with anaesthetic agents.

4.5 Filling Adapter

For safe filling of the Sedaconda Syringe the correct filling adapter must be used.

4.6 Gas scavenging system

Sedana Medical recommends scavenging the exhaust gases from the ventilator and the gas monitor.

Passive Gas Scavenging

There is a passive scavenging system available from Sedana Medical called FlurAbsorb and FlurAbsorb-S which is used in conjunction with an accessory kit.

Active Gas Scavenging

Active Gas Scavenging can be used if installed or a central vacuum source can be used together with a pressure equalization system, which can be provided by the manufacturer of ventilators.


5. SYSTEM ASSEMBLY

5.1 Gas scavenging system

- Connect the exhaust from the ventilator and gas monitor to the gas scavenging system.

5.2 Filling the Sedaconda Syringe

- Open the bottle and screw the filling adapter on to the bottle.
- Unscrew the red cap from the red top of the syringe.
- Aspirate air in the syringe and connect the syringe tightly to the filling adapter.
- Turn the bottle upside down. Fill the syringe by slowly moving the syringe plunger back and forth.
- Turn the bottle and syringe horizontal upright and unscrew the syringe from the filling adapter.
- Remove any air from the syringe and close the syringe with the red cap.
- Note on the syringe label which anaesthetic agent is used and date of filling.

5.2. Set-up (Fig.1)

- Ensure that the syringe pump is at or below the level of the Sedaconda ACD.
- Remove the red protective cap from the Sedaconda ACD.
- Connect the gas sampling line to the Sedaconda ACD gas sampling port.
- Connect the other end of the gas sampling line to the gas analyser port.
- Insert the Sedaconda ACD in the breathing circuit between the ET tube and the Y-piece. Ensure the Sedaconda ACD has the black side up and is sloped down towards the patient.
- Open the red cap on the Sedaconda Syringe and connect the Sedaconda ACD anaesthetic agent line to the syringe.
- Place the filled Sedaconda Syringe in the syringe pump.

6. OPERATING

6.1 Priming the agent line

- Administer a bolus of 1,2 ml when initially connecting Sedaconda ACD. If necessary, give an extra bolus of 0.3 ml. Never give manual bolus.
- Set the clinical dosage
- Start the syringe pump (check point 6.2 below)

**6.2 Dosing the Anaesthetic agent**

All dosing is individual and guided by experienced clinical evaluation and hemodynamics. There is a higher patient uptake of the volatile during the first 10-30 minutes (Induction Phase) of administration and therefore corrections of the pump rate need to be made according to the clinical needs of the patient. Isoflurane is approximately twice as potent as sevoflurane.

The following rates are typical for the initial syringe pump rate of isoflurane and sevoflurane:  
 – isoflurane: 3 ml/h – sevoflurane: 5 ml/h

The syringe pump rate necessary to reach a certain patient concentration depends on the minute volume.

Volatile Agent	Expected Pump Rates	Resulting Fet Values
isoflurane	2 – 7 ml/h	0.2 – 0.7%
sevoflurane	4 – 10 ml/h	0.5 – 1.4%

If a rapid increase of the concentration is deemed necessary, a bolus of 0.3 ml may be given. At higher tidal volumes and / or higher respiratory rates, the Sedaconda ACD is less efficient. Therefore relatively more anaesthetic, and thus a higher pump rate is needed to keep the concentration stable.

**6.4 Ending the Therapy**

**Immediate Cessation**

1. Stop the syringe pump.
2. Disconnect the agent supply line from the Sedaconda Syringe.
3. Seal the syringe with the syringe closure.
4. Disconnect the gas monitor from the Sedaconda ACD. Close the gas monitor port with the gas sampling port closure.
5. Remove the Sedaconda ACD from the patient. Disconnect from the Y-piece first.
6. The concentration will decrease rapidly.
7. Consider replacing the Sedaconda ACD with a bacterial/viral filter with heat and moisture exchanger.
8. Dispose of the Sedaconda ACD according to hospital protocol.

**Short Weaning Process**

1. Stop the syringe pump and leave the Sedaconda ACD in place.
2. The concentration will gradually decrease.
3. To accelerate the final steps in the weaning process follow above steps (2-8) under 'Immediate Cessation'.

**6.5 Changing the Sedaconda ACD**

- Prepare a new Sedaconda ACD, and a new filled syringe if needed (as per 5.1).
- Stop the syringe pump.
- Disconnect the anaesthetic agent line from the syringe and close the syringe with the red cap.
- Disconnect the gas sampling line from the old Sedaconda ACD and close the gas sampling port.
- Connect the gas sampling line to the new Sedaconda ACD
- Disconnect the old Sedaconda ACD from the Y-piece first, and then from the ET tube.
- Insert the new Sedaconda ACD. Ensure the Sedaconda ACD has the black side up and is sloped down towards the patient.
- Connect the agent line to the syringe in the syringe pump.
- Secure all connections.
- Prime the agent line with 1.2 ml. Prime with the syringe pump. Never prime manually.

**6.6 Changing the Sedaconda Syringe**

- Stop the syringe pump.
- Remove the empty syringe from the syringe pump.
- Disconnect the anaesthetic agent line from the syringe and close the syringe with the red cap.
- Unscrew the red cap and connect the Sedaconda ACD agent line to the new syringe.
- Place the full syringe in the syringe pump.
- Start the syringe pump with the same rate as before.
- Dispose of the old syringe according to hospital protocols.

**7. CONNECTING A NEBULIZER TO THE SEDACONDA ACD SYSTEM**

It is possible to use a jet or ultrasonic nebulizer with the Sedaconda ACD system. The nebulizer should be connected between the patient ET tube and the Sedaconda ACD. Ultrasonic nebulizers are preferable as they do not add extra airflow. If a jet nebulizer is connected it may be necessary to increase the syringe pump rate, to compensate for the extra flow from the nebulizer. When connecting a nebulizer set the ventilator on stand-by or hold an expiratory pause on the ventilator.

**Repeated nebulizations may increase the flow resistance of the Sedaconda ACD. Pay attention to signs of occlusions. Sticky nebulized drugs (e.g. acetylcysteine, colistin and amphotericin B) may increase the resistance of the filter and require early replacement of the Sedaconda ACD.**

**NOTE!** Always consider the increased dead space when connecting extra items.

**8. SUCTIONING**

- Using a closed suction system or using a swivel connector with suction port is preferable.
- Hold pause on the ventilator if disconnecting Sedaconda ACD from the ET-tube during the procedure. When disconnecting, remove the Sedaconda ACD from the Y-piece first and when attaching, attach the Sedaconda ACD to the ET-tube first.

**Please be aware that devices or components used in the patient breathing circuit which are composed of polycarbonate or polystyrene may become degraded or undergo stress cracking in the presence of isoflurane or sevoflurane anaesthetic gases.**

**9. DISPOSAL**

Dispose of the Sedaconda ACD and the sealed Sedaconda Syringe according to hospital protocols.

**10. TECHNICAL INFORMATION**

TECHNICAL SPECIFICATION	Sedaconda ACD-S	Sedaconda ACD-L
Anaesthetic Agents	Only use room temperature sevoflurane (18°-25°C) and isoflurane (18°-30°C)	
Syringe	Only use the Sedaconda Syringe	
Stability of filled syringes	5 days	
Tidal volume working range	>200 ml	>350 ml
Sedaconda ACD dead space	Approx. 50 ml	Approx. 100 ml
Moisture loss	5 mg/l (@ 0.5L X 15 bpm) 6 mg/l (@ 0.75L x15 bpm)	5 mg/l (@ 0.75L X 12 bpm) 7 mg/l (@ 1.0L x 10 bpm)
Compliance @ 70cm H2O	1 ml/kPa	2 ml/kPa
Gas leakage*	< 0.01 L/min	
Bacterial Filtration Efficiency	99.867%	
Viral Filtration Efficiency	99.76%	
Weight	40 g	52 g
Agent Line Length	2.2 m	
Connectors (According to ISO 5356)	Patient Side Connector: 15F/22M Ventilator Side Connector: 15M	
Gas Sampling Port	Female Luer Taper Connector	

\*During the device entire usage period

**11. SEDACONDA ACD – VENTILATOR/ALTERNATIVE PLACEMENT**

The Sedaconda ACD is inserted on the inspiratory limb of the breathing circuit, by the ventilator, positioned with a downward slope.

This alternative placement is intended for, and enables the delivery of inhaled anaesthetics to patients with small tidal volumes (30-200 ml) or other situations where dead space/carbon dioxide accumulation is an issue.

Always use the large FlurAbsorb and replace it after 5 syringes (each 50 ml) when using active humidification.

If active humidification is being used place the heater below the Sedaconda ACD, to avoid accumulation of condensate.

**Do Not** use the Sedaconda ACD gas sampling port, since it is not giving correct gas-measurements in the ventilator/alternative placement. Measuring of gas concentrations should be done at the Y-piece.


Only use active humidifiers that are compatible with volatile agents. Replacement of the Sedaconda ACD and/or Sedaconda Syringe should be done swiftly, to minimise time without drug delivery, since there is no gas reflection in this ACD placement. Always take the age and size of the patient into account when giving bolus.

In contrast to the standard placement of the Sedaconda ACD (between the ventilator/Y-piece and the patient), the ventilator/alternative placement only utilizes the evaporator function of the Sedaconda ACD and hence there is no reflection of inhaled anaesthetic. Therefore, higher anaesthetic pump rates can be expected, as there is no reflection of anaesthetic agent.

There is no reflection of heat and moisture either, therefore other means of humidification is necessary.

**Materials needed for the assembly (fig 2)**

Fig. 2

**In case of a serious incident please contact the competent authority and the legal manufacturer Sedana Medical Ltd.**

For further information regarding policies or procedures relating to the Sedaconda ACD the user should contact Sedana Medical AB.

Sedana Medical Ltd., Unit 2A The Village Centre,  
Two Mile House, Naas, Co. Kildare, W91 PWH5, Ireland

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