Čištění

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English

Notification of Adverse Events

As a health care provider, you may report the occurrence of certain events to SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD., and possibly to the competent authority of the Member state in which the user and / or patient is established.

These events, include device-related death and serious injury or illness. In addition, as part of our Quality Assurance Program, SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. requests to be notified of device failures or malfunctions. This information is required to ensure that SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. provides only the highest quality products.

Intended use

Intended purpose

The Artema Dryline Water traps are intended for docking with the Artema Dryline Receptacle and connecting to an Artema Dryline Sampling Line in a side-stream gas monitoring system. The Water traps filtrates the gas sample from humidity, secretions, bacterial contamination and dust. There is one Dryline Water trap dedicated for Adult/Pediatric use. For use only by health-care professionals trained in anesthesia or intensive care.

Medical conditions

For use in a side-stream gas monitoring system.

Intended users

For use only by health-care professionals trained in anesthesia or intensive care.

Intended patient population

The products can be used in adult and pediatric patients.

Contra-indications

There are no absolute contraindications for this product.

Side effects

N/A

According to the conclusion of clinical evaluation and residual risk evaluation, for the intended patients, there is no known side effects that can occur during or after the use of the medical device. And there is no need for the operator to make extra preparations. Thus, no residual risk associated with using the medical device should be disclosed due to the risk management report.

Instructions before use

- To install the water trap, align it with the receptacle and push gently into place (Figure 1). Make sure that the locking mechanism is fully engaged by pulling the water trap, which should be firmly seated.
- Connect an Artema Technology DRYLINE™ Sampling Line, Adult (color-less tube) to the gas sampling port.
- Perform a pre-use test (refer to the host instrument Instruction For Use).

Warning

- No modifications of this equipment is allowed.
- Connect only Artema Technology adult gas sampling lines to the water trap. Other tubes, e.g. IV-lines, may be possible to attach but doing so could result in patient harm.
- Do not use with Artema Technology DRYLINE™ Sampling Line, Neonate (blue tube).
- Handle the contents of the water trap as you would handle any body fluid.

Emptying/replacing the water trap

- Maximum emptying interval during normal use (37°C, 100% RH) is 17 hours @200ml/min and 26 hours @120 ml/min.
- To remove the water trap, push the release button on the receptacle and pull out the water trap. (Figure 2).
- To empty the water trap, twist and pull the container relative to the filter housing. Empty the container (Figure 3 & 4) and re-install the water trap (Figure 1).
- Replace the complete water trap every month or more often if the host instrument indicates this.

Cleaning

- No parts of the water trap are intended to be cleaned.

Warning

- Never allow detergent to enter the filter housing.
- Never force air through the water trap.

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