


Gas Sampling Line

Product description

The Gas Sampling Line is for use with Sedaconda® ACD, please refer to the Sedaconda ACD Instructions for Use, and Bistos BT-740.



Features & Benefits

- Approved with isoflurane, sevoflurane and Sedaconda ACD
- Unique appearance and length
- Remove moisture and humidity
- Simplify the set-up
- Reduce water in the water trap
- Easy to differentiate from other medical lines

- **Device Classification:** Class I
- **Use:** Single Usage 72 h
- **Shelf Life:** 3 Years from Date of Manufacturing

- **Reference Number:**
- **Disposal:** According to hospital protocol
- **Packaging Unit:** 1 unit contains 12 Gas sampling lines
- **Length:** 2.5m/8ft

Distributor

Sedana Medical Ltd.
Unit 2A Village Centre
Two Mile House Naas,
Kildare, Ireland

Legal Manufacturer



For questions, please contact Sedana Medical at medinfo@sedanamedical.com

For product feedback and complaints, please contact: safetyandcomplaints@sedanamedical.com

Before use, please refer to the Instructions for Use (IFU) for each product

XXXXXXX